

New – StatSafe Drawer Sub-dividing



New – StatSafe Drawer Sub-dividing

The new StatSafe Drawer Sub-dividing features allow you to sub-divide open matrix drawers with up to 10 subsections. This feature can be used with any open matrix drawer or bin drawer location.

Dividers can be used to create up to 10 sub-sections within each row in the drawers. You can see our drawer organization accessories here:

https://phoenixltc.com/product-category/accessories/drawer-organizers/



To set up this new feature, follow these steps:

- Your StatSafe will need to be on an updated version of software to see these features. If these do not show as available in the StatSafe Portal, at the StatSafe login screen, select "Options" and "Update Program" to update the StatSafe to the latest software version.
- 2. In the StatSafe Portal, select the StatSafe you would like to apply these features to and select "Edit Settings". Then, under the "Machine Config" tab, select the drawer that you would like to add the sub-divided sections to.

- **3.** Under "Config Options" choose "Allow Row Sub-dividing".
- After the previous steps are complete, whenever you create inventory in drawer locations with sub-dividing turned on, you will be prompted to choose from slots 0-9.

Note: slot 0 is not an assigned slot and just adds the items to that drawer and row.

Select a Location		
ACTIVE		
Item Location Drawer 3 section C		
(Space is available here)		
Slot		
Slot 3		*
Items in this location		
	Search	Q
ID		
No records to display		

5. The labels will now show the drawer, row, and slot on them.

Config Options			
	Dispense only from this location if Rx available		
	Dispense from this location if Rx available but allow override		
	Allow all dispensing from this location even if StatSafe is set for enforce orders		
	Allow dispensing from this location from Rx Orders even if StatSafe has been offline > 8 hours		
	Disable witness when accepting		
	Perform shift audits on this location		
	Allow patient assigned items		
	Patient assigned items only		
	Don't send patient assigned items to interface		
	Disable blind count on this location		
\checkmark	Allow Row Subdividing		



6. When a nurse accepts stock, it will now show the drawer, row, and slot to place the item in.

	3B-720 Slot 3 Skittles Drawer 3	ACCEPT STOCK Please follow the steps below to accept this inventory item into the system
STEP 1	Pin (Found on back or bottom o	
STEP 2	Quantity (Confirm Qty)	
STEP 3	Load Item (Place item into location)	
STEP 4	FINISH	Cancel 🙁

7. When a nurse pulls the medication, it will also show on the screen which slot to choose from:

