

STATSAFE

Automated Medication Management

By  PHOENIX



StatSafe – New Two-Factor Authentication



PhoenixLTC.com | 855.MED.CART (633.2278)

StatSafe – New Two-Factor Authentication


StatSafe now allows for two-factor authentication for logging into the StatSafe medication management systems.

There are three methods of user access to log in to StatSafe.

**Biometrics are not available on all StatSafes. Please reach out if you would like more information on our latest StatSafes with Biometric login.

1. Username and login
2. Biometric login
3. Card verification using a scanned barcode registered to the user

To set up two-factor authentication, select the StatSafe you would like to apply the settings to. Under the cart settings and S2 settings, select "Use 2 Factor Authentication on Login" under "Additional Options". Then, under "Security Options" choose the methods to be used by the facility users.



The screenshot displays the StatSafe configuration interface. It is divided into three main sections:

- Adjust StatSafe Time to:** A row of radio buttons for selecting a time interval: -1 Hr, 0 Hr (selected), 1 Hr, 2 Hrs, and 3 Hrs.
- StatSafe Password Timeout (Days):** A numeric input field with the value '0' and up/down arrow controls.
- Security Options:** A list of checkboxes where 'Enable Card Verification' and 'Enable Biometric ID' are both checked.
- Additional Options:** A list of checkboxes on the right side of the interface. The checkbox for 'Use 2 Factor Authentication on Login' is checked, while all other checkboxes are unchecked.

When the users log in, they will be prompted to enter their username and password, then to provide their secondary login method after that.

When these features are turned on, it must apply to all users.